Weropress® / Weropress® LT

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Weropress® / Weropress® LT

Cold-curing polymer

(en) Instructions for use, please read carefully!

Intended use

Cold-curing, colour-stable, methylmethacrylate-based resin for fabricating partial and total prostheses using injection, tamping-pressing and casting processes.

Indication

- Total or partial tooth loss
- Abnormal jaw function

Contraindication

- Hypersensitivity to one of the components

Composition

| Weropress® Polymer | Weropress® Monomer | Weropress® LT Monomer |
|---------------------------------|---------------------------------|---------------------------------|
| PMMA (polymer und copolymer) | MMA (methylmethacrylate) | MMA (methylmethacrylate) |
| barbituric acid catalyst system | dimethacrylate | dimethacrylate |
| organic colourants | barbituric acid catalyst system | ethermethacrylate |
| inorganic colourants | | barbituric acid catalyst system |

Injection process

- Proceed as specified by the device manufacturer
- The temperature of the plaster in the flask should be 35 - 40 °C / 95 - 104 °F to promote polymerisation.

Press process

- Flask technique (injection and tamping-pressing processes)
- The temperature of the plaster in the flask should be 35 - 40 °C / 95 - 104 °F to promote polymerisation.
- Complete the press process promptly under the hydraulic press within 1 minute.

Casting process

- Investment with plaster, silicone or hydrocolloid (gel) depending on the process.
- Polymerisation is conducted at 2 to 2.5 bar pressure in water at 45 °C / 113 °F after appropriate dwell time.


Dosage and application

- See table for mixing ratios and processing times depending on the processing method of the acrylic.
- Use class IV plaster for making casts and investment where possible.
- Always wet casts sufficiently for 5 to 15 minutes depending on the degree of dryness.
- Use standard alginate isolation to separate plaster from acrylic.
- Do not compress the powder when using dosage aids.

Times and quantities

| | Injection process | Tamping-pressing process | Casting process |
|-----------------------------------------|----------------------------|----------------------------------------------------------------|-------------------------------------------|
| Mixing ratio* | 30 g : 15 ml | 30 g : 15 ml | 30 g : 21 ml |
| Polymerisation temperature | Residual heat in the flask | Residual heat in the flask | 45 °C / 113 °F |
| Polymerisation time | 15 min | 15 min | 15 min ¹ / 30 min ² |
| Pressure required during polymerisation | By device type | Pressure in the flask must be maintained during polymerisation | 2 - 2,5 bar |

Weropress®

Weropress® LT

Light grey box: Mixing time

Dark grey box: Dwell time

Black box: Processing time (pourable consistency)
(Start and end of injection time/tamping time)

¹ in matrix technique, plaster, putty

² in flask duplication technique, duplication gel,

*see also dosage aid **REF 1090194**

⁴ Based on a room temperature of 21 °C / 69,8 °F

⁵ Based on a room temperature of 30 °C / 86 °F

Trimming

To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

Bond

- Roughen the basal surfaces of the acrylic teeth before bonding.
- Ensure that no impurities such as wax or isolating agent are present.
- Mechanical retention is recommended.
- Adequate mechanical retention is essential for ceramic teeth.

Restoration

- Products of Weropress® / LT can also be restored and supplemented at any time with all standard MMA-based cold-curing polymers.
- Following roughening of the surface (sanding / blasting / grinding), further conditioning is recommended, preferably with artConnect PLUS or a similar MMA-based bonding agent.

Disposal

- The fully polymerized components of Weropress® / LT are disposed of as residual waste.
- Disposal of nonpolymerized product residues according to legal regulations.
- Do not allow entry into the sewage system or waterways.

Notes

- Close all powder and liquid containers tightly after use.
- Store in a dry and cool area away from direct sunlight.
- Changes to the construction through use, e.g. due to fissures, cracks, clouding or significant signs of abrasion, necessitate a repair or refabrication. Monitoring should be carried out with regular assessment of the denture at the follow-up presentation in the dental surgery.
- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.
- Processing the acrylic: To reduce the residual monomer content, the finished (sanded) but not polished polymerised workpiece should be stored in water at room temperature for about 48 hours before delivery.
- The residual monomer content after polymerisation of Weropress® / LT is < 4,5% (EN ISO 20795-1).

Warnings

Processing frameworks may generate dust which can lead to mechanical irritation of the eyes and airways. Therefore, always ensure that the extraction system on your milling machine and at the workplace is working properly for individual post-processing and that you are using personal protective equipment.

Danger and safety instruction Weropress® & Weropress® LT Monomer

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P233 Keep container tightly closed.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P262 Do not get in eyes, on skin, or on clothing.

P210 Keep away from heat. No Smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

Danger and safety instructions Weropress® Polymer

EUH208 Contains dibenzoyl peroxide; benzoyl peroxide, methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate. May produce an allergic reaction.

Classification

Weropress® / LT conforms to EN ISO 20795-1 Type 2 Class 1.

The product properties are based on compliance with and observation the instructions for use.

For detailed information please see the applicable Material Safety Data Sheet, available upon request or as download from www.merz-dental.de.

Date of information 2020-06

Weropress®

Cold-curing polymer

USA **CAN** **Instructions for use, please read carefully!**

Caution: Federal law restricts this device to sale by or on the order of a licensed dental practitioner.

Cold-curing, colour-stable, methylmethacrylate-based resin for fabricating partial and total prostheses using injection, tamping-pressing and casting processes.

Composition

Weropress® powder contains

- PMMA (polymer and copolymer)
- barbituric acid catalyst system
- organic colorants
- inorganic pigments

Weropress® liquid contains

- MMA (methylmethacrylate)
- dimethacrylate
- barbituric acid catalyst system

Dosage and application

- See table for mixing ratios and processing times depending on the processing method of the acrylic.
- Always wet casts sufficiently for 5 to 15 minutes depending on the degree of dryness.
- Use standard alginate isolation to separate plaster from acrylic.
- Do not compress the powder when using dosage aids.

Injection process

- Proceed as specified by the device manufacturer
- Use class IV plaster for making casts and investment where possible.
- The temperature of the plaster in the flask should be 35 - 40 °C / 95 - 104 °F to promote polymerisation.

Press process

- Use class IV plaster for making casts and investment where possible.
- The temperature of the plaster in the flask should be 35 - 40 °C / 95 - 104 °F to promote polymerisation.
- Complete the press process promptly under the hydraulic press within 1 minute.


Casting process

- Use class IV plaster for making casts where possible.
- Investment with plaster, silicone or hydrocolloid (gel) depending on the process.
- Polymerisation is conducted at 2 to 2.5 bar pressure in water at 45 °C / 113 °F after appropriate dwell time.

Times and quantities

| | Injection process | Tamping-pressing process | Casting process |
|-----------------------------------------|----------------------------|----------------------------------------------------------------|-------------------------------------------|
| Mixing ratio* | 30 g : 15 ml | 30 g : 15 ml | 30 g : 21 ml |
| Polymerisation temperature | Residual heat in the flask | Residual heat in the flask | 45 °C / 113 °F |
| Polymerisation time | 15 min | 15 min | 15 min ¹ / 30 min ² |
| Pressure required during polymerisation | By device type | Pressure in the flask must be maintained during polymerisation | 2 - 2,5 bar |

Processing the acrylic

To reduce the residual monomer content, the finished (sanded) but not polished polymerised workpiece should be stored in water at room temperature for about 48 hours before delivery.

Trimming

To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

Bond

- Roughen the basal surfaces of the acrylic teeth before bonding.
- Ensure that no impurities such as wax or isolating agent are present.
- Mechanical retention is recommended.
- Adequate mechanical retention is essential for ceramic teeth.

Restoration

Products of Weropress® can be restored and supplemented at any time with all standard MMA-based cold-curing polymers, such as Weropress®.

Notes

- Close all powder and liquid containers tightly after use.
- Store in a dry and cool area away from direct sunlight.
- Take precautionary measures against static discharges.
- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.
- The residual monomer content after polymerisation of Weropress® is < 4,5% (EN ISO 20795-1).

Disposal information

Disposal of contents/container in accordance with local/regional/national/international regulations.

Danger and safety instruction Weropress® Monomer

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P233 Keep container tightly closed.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P262 Do not get in eyes, on skin, or on clothing.

P210 Keep away from heat. No Smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

Danger and safety instructions Weropress® Polymer

EUH208 Contains dibenzoyl peroxide; benzoyl peroxide, methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate. May produce an allergic reaction.

Classification

Weropress® conforms to EN ISO 20795-1 Type 2 Class 1.

The product properties are based on compliance with the instructions for use.

For detailed information see the applicable Material Safety Data Sheet.

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Merz Dental GmbH